**Department of Civil Engineering,
Maulana Azad National Institute of Technology,
Bhopal**

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Date/ दिनांक.....

Objective: Verification of Bernoulli's theorem.

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Apparatus required:

1. Apparatus for verification of Bernoulli's theorem.
2. A collecting tank
3. A stop watch
4. A meter scale

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Theory:

The Bernoulli's theorem states that in a steady continuous flow of an ideal incompressible fluid, the sum of the velocity head, the pressure head and the potential head is constant i.e. Each term of equation describes energy in Kg-m per Kg and has the dimensions of length. The terms $(v^2/2g)$ is known as the velocity head, (p/γ) is known as the pressure head and Z as the potential head. Applying Bernoulli's equation between any two points, one gets,

$$\frac{v_1^2}{2g} + \frac{p_1}{\gamma} + z_1 = \frac{v_2^2}{2g} + \frac{p_2}{\gamma} + z_2$$

Hence, if V_2 is greater than V_1 , $(Z_1 + (p_1/\gamma))$ must be greater than $(Z_2 + (p_2/\gamma))$. In other words, any increase in velocity will cause a corresponding reduction in the $(Z + (p/\gamma))$. Similarly, a reduction in velocity will cause an increase in $(Z + (p/\gamma))$. In case the flow is in a horizontal plane or when the fluid is a gas, the term Z drops out and hence changes in velocity cause corresponding changes in the pressure.

In real fluids, which have viscosity, some energy is converted into heat due to the viscous shear and consequently there is a certain amount of energy loss. Hence, for real fluids equation can be written as follows:

$$\frac{v_1^2}{2g} + \frac{p_1}{\gamma} + z_1 = \frac{v_2^2}{2g} + \frac{p_2}{\gamma} + z_2 + h_L$$

In which the subscript 1 and 2 refer to the two sections under consideration and h_L is the energy loss per unit weight of fluid between the two sections.

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Procedure:

1. Open the inlet valve to reservoir so that small quantity of water starts flowing into the reservoir
2. Regulate the control valve in pipeline so that the water level in the reservoir remains constant for an opening of inlet valve.
3. Find the discharge flowing by volumetric measurement in the collecting tank by noting the rise of water level in the collecting tank for a known interval of time.
4. Record the levels in all the piezometers and the cross-section of the channel at each piezometer.
5. Change the opening of the inlet valve and repeat the steps 2, 3 and 4.
6. Take 2-3 sets of readings.

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Calculations:

1. Calculate the discharge from volumetric measurements.
2. Find the velocity of flow at all the piezometric points by dividing the discharge by the corresponding area of channel and then calculate the velocity head at these points.
3. Find the sum of all the 'energies at each point for a given flow. Repeat the calculations for all sets of readings.
4. Tabulate the results.

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Presentation of results:

1. On a simple graph, mark the position of tube on the abscissa and $(p/\gamma + Z)$ on the ordinate. Join all ordinates to give hydraulic gradient line (HGL).
2. Mark velocity head, $V^2/2g$ over $(p/\gamma + Z)$. Join the points to make total energy line (TEL)
3. Add loss of head to TEL and draw a horizontal line showing theoretical TEL.

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Questions

1. What are the assumptions made in deriving the Bernoulli's equation?
2. What are the applications of Bernoulli's equations?
3. What are steady and unsteady flows?
4. What additional terms will be introduced in the present form of Bernoulli's theorem if the flow is unsteady and incompressible?

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Date / दिनांक: _____

Objective: To verify Impulse-Momentum principle for flow over vane.

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Apparatus required

1. An experimental set up consisting of nozzle issuing upward vertical jet of water against a vane under test fixed on a vertical rod with a lever arm to carry load
2. A collecting tank
3. A stop watch
4. Weights or loads

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Theory:

Many problems cannot be easily solved by consideration of energy principle (Bernoulli's equation) alone due to complicated geometry of flow in some situations. Examples of such flows are in pumps, turbines, vanes, pipe bends & enlargements etc. In such cases impulse momentum principle can be conveniently used. It states that the impulse of the force or product of the force and the time increment during which it acts is equal to the change in the momentum of the body i.e. product of the mass of the body on which the force acts and the velocity of the body.

If F_x and F_y are the components of the force F along the two reference coordinate directions, V_x and V_y the corresponding velocity components and Q the discharge.

$$F_x = \frac{Q \cdot \gamma}{g} (V_{2x} - V_{1x}), \quad F_y = \frac{Q \cdot \gamma}{g} (V_{2y} - V_{1y}) \quad \text{and} \quad F = \sqrt{F_x^2 + F_y^2}$$

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Procedure:

1. Fix a vane to the vertical rod. Note the weight of the vane, the rod and collar etc. fixed to it. Alternately, balance the weight of the vane, rod and the lever by sliding a balancing weight beyond the fulcrum of the lever.
2. Open the valve controlling the jet and adjust discharge such that the force of the jet is just equal to the total weight (W) acting down on the vane and the vane and the lever are balanced. Give a slight motion to the lever to minimize friction.
3. Note the discharge in the collecting tank.
4. Repeat by increasing the jet velocity and weight on vane.
5. Take different sets of readings.
6. Repeat using different vanes. Note the deflection angle of each vane.

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Calculations:

Calculate the rate of flow through the nozzle, $Q = (Axh)/t$ and then velocity of jet, $V = Q/a$

$F_{\text{theoretical}} = \rho QV(1+\cos\theta)$, $F_{\text{actual}} = w L_2/L_1$ g, $\% \text{Error} = (F_{\text{theoretical}} - F_{\text{actual}})/F_{\text{actual}}$

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Sample Calculation / नमूना गणना:

Presentation of results:

Plot F as abscissa Vs F_a as ordinate. One should get a straight line at 45° to the horizontal.

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Date / दिनांक: _____

Objective:

To calibrate the venturi meter and study the effect of Reynolds number on venturi meter coefficient, C_d .

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Apparatus required

1. Venturi meter fitted in a pipe line
2. Differential monometer connected to the upstream enlarged end and throat of the venturi meter
3. A collecting tank
4. A stop watch
5. A thermometer

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Theory:

Venturi meter is a rate meter used for measuring flow through closed conduits, circular in cross-section and is especially suited for measuring large flows with least head loss. It is considered to be far more superior to other meters like orifice meter and flow nozzle of the class though more costly. It consists essentially of a convergent divergent streamlined tube. The streamlined boundary contraction and expansion cause a gradual pressure rise in the divergent section. To avoid separation of flow and the consequent loss of energy due to formation of eddies, the angle of convergence is limited to a maximum of 20° and the angle of divergence to a maximum of 7° . The present-day venturi meters have a short cylindrical throat section of diameter from

1/4 to 3/4 (usually 1/2) of the pipe diameter. For accurate results, the venturi meters should be preceded by a straight pipe of at least 30 pipe diameters.

To find the discharge through the venturi meter, Bernoulli's equation may be applied to pipe section (1) and the throat section (2). Then

$$\frac{v_1^2}{2g} + \frac{p_1}{\gamma} + Z_1 = \frac{v_2^2}{2g} + \frac{p_2}{\gamma} + Z_2$$

where Z_1 and Z_2 refer to the distance of the centre of the sections of the venturi meter above an arbitrary datum and V_1 and V_2 and p_1 and p_2 refer to the uniform velocities and pressures at sections (1) and (2) respectively. Also, by continuity equation

$$A_1 V_1 = A_2 V_2$$

Solving these two equations:

$$V = \frac{1}{\sqrt{\left\{1 - \left(\frac{d_2}{d_1}\right)^4\right\}}} \times \sqrt{2g \left(\frac{p_1 - p_2}{\gamma}\right)}$$

Multiplying by cross-sectional area A_2 of throat, theoretical discharge,

$$Q_{th} = A_2 V_2 = \frac{A_2}{\sqrt{1 - \left(\frac{d_2}{d_1}\right)^4}} \times \sqrt{2g \left(\frac{p_1 - p_2}{\gamma}\right)}$$

$$Q_{th} = \frac{a_1 a_2 \sqrt{2gh}}{\sqrt{a_1^2 - a_2^2}}$$

If a differential manometer is now connected to points (1) and (2) and if h is the gauge reading, then by manometry principles

$$(p_1 - p_2) / \gamma = h \left\{ \frac{s_1}{s} - 1 \right\}$$

Where s_1 and s are specific gravities of gauge liquid and fluid flowing.

Further there will be loss of head between the sections (1) and (2) and this tends to increase the observed value of h . To take this into account, a coefficient C_d called the coefficient of discharge is introduced thus modifying the discharge equations of the venturi meter to C_d is usually obtained by calibration of the meter.

$$Q_{act.} = \frac{C_d A_2}{\sqrt{1 - \left(\frac{d_2}{d_1}\right)^4}} \times \sqrt{2g \left(\frac{p_1 - p_2}{\gamma}\right)}$$

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Procedure:

1. Note the dimensions of the cross section of collecting tank and sections of the venturi meter;
2. Keep the drain cock of the collection tank open and open valve on the pipe line partly to allow some flow through the venturi meter. Open the cocks of the venturi meter gradually and simultaneously thus connecting up the manometer. Open the flushing cocks of the manometer gradually and simultaneously to remove entrapped air bubbles in the manometer leads and after flushing close the flushing cocks.
3. Note down the reading of the monometer and temperature of water.
4. Close the drain cock of the collecting tank. The water level in the tank rises. When the level has reached a convenient full scale reading on the scale attached to the gauge glass of the tank, start the stop watch and note the time taken for the known rise of water level in the collecting tank. This completes one set of readings.
5. Vary the discharge through the Venturi meter by manipulating the valve on the pipe line Allow the condition of flow to become steady. Note the manometer readings and the time for the water level in the collecting tank to rise by a known height after

lowering the water level in the tank to a convenient level, if necessary, by opening the drain cock.

6. Repeat the experiment with various discharges to obtain at least ten sets of readings.
7. Close the valve on the upstream side of the venturi meter completely, thus cutting off the flow. Close the cocks of the venturi meter.

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Calculations:

1. Calculate the area of cross-section of the collecting tank and the area A2 of the throat section of the venturimeter.
2. From the area of cross-section of the collecting tank and the observed time for water level in the tank to rise by a known height, calculate the discharge, Q for each of readings.
3. From equation (6) calculate the values of the coefficient of discharge, Cd for each set of observations Si and S being known or assumed.
4. Calculate Reynolds Number for each discharge as:

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$$Re = \frac{V_2 \times d_2}{\nu}$$

Representation of Results:

Draw the curves (1) manometric head, h Vs discharge, Q on the double log graph paper (2) C_d with Reynolds number on semi-log graph paper.

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Questions:

1. What are the advantages and dis-advantages of venturi meter as compared to other similar meters? When do you use a venturi meter?
2. What happens if the convergence and divergence angles of a venturi meter are increased?

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Date / दिनांक: _____

Objective:

To calibrate and orifice meter and determine the variation of orifice meter co-efficient with Reynolds number

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Apparatus Required:

1. Orifice meter filled in a pipe line
2. Differential manometer connected to upstream section and at vena contracta
3. A collecting tank
4. A stop watch
5. Thermometer

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Theory:

An orifice meter is another device used for the measurement of flow rate in a pipe line. It consists of thin metal plate with an opening (generally circular concentric with the pipe axis). The plate is held between two flanges in the pipe line. The flow characteristics of an orifice meter differ from those of a venturi meter. The stream lines approaching contracts and diverge suddenly before and after orifice meter. The minimum cross- section of the jet occurs not within the orifice but downstream from the orifice edge (approximately at a distance of half the orifice diameter). The section of minimum area (a_c) is called vena-contracta. The ratio $C_c = a_c/a_0$ is called coefficient of contraction Discharge co-efficient of an orifice meter varies with the type of orifice the pipe size, ratio of orifice diameter to pipe diameter., Reynolds number and location of pressure connections. The device developed based on the concept of Bernoulli's theorem. As per the Bernoulli's theorem, for constant potential energy and flow rate, increase in velocity will decrease the pressure at any point and this concept can be used for discharge measurement in pipes. Orifice meter consists of a. orifice plate as shown in above

figure. The amount of pressure changes between throat and inlet depending on the flow rate. As flow area is suddenly reduced hence due to losses actual discharge will be less than theoretical discharge.

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Pressure head in terms of water (h) $h = 12.6 h_m$

Coefficient of Discharge $C_d = Q_a / Q_t$ Where, a_1 – cross sectional area of inlet pipe a_2 – cross sectional area of orifice A – area of measuring tank y – rise of water level in measuring tank during the time t h_m – difference of mercury columns in differential U tube manometer

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Procedure:

1. Open the manometer cocks and ensure that there are no air bubbles present in connecting tubing
2. Open the outlet valve so that water blows through the pipe line and past orifice meter wait for some time so that the flow is stabilized,
3. Note the manometer reading in the two limbs of the manometer
4. Measure the discharge by collecting water to a certain depth in the collecting tank for some time interval 't'
5. Change the valve opening and repeat steps 2-4 for 6 to 8 sets of reading

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Calculations:

Following calculations are to be carried out

1. Actual discharge, $Q_a = (Axh)/t$.
2. Theoretical discharge Q_t .
3. Coefficient of Discharge C_d
4. Finally, calculate Reynolds number, $Re = \frac{V_1 \times d_1}{\nu}$ corresponding to different discharges.

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Presentation of Results:

Plot the following graphs:

1. Plot $Q_{act} V_s h$ on a double log graph paper and fit in straight line for the plotted points.
2. Plot $C V_s R_e$ on a semi- log graph.

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Questions:

1. What is use of an orifice meter? Compare it with a venturi meter?
2. Is there any factor to limit the minimum value of d/D that is workable?
3. Why is vena-contracta formed?

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Date / दिनांक: _____

Objective:

To calibrate the nozzle meter and find its coefficient of discharge.

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Apparatus Required:

1. Nozzle meter fitted in pipe line.
2. Differential manometer connected to the upstream of the nozzle and a point just downstream of the nozzle.
3. A collecting tank.
4. A stop watch.
5. A thermometer.

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Theory:

Nozzle meter is a flow meter used for measuring flow through closed conduits. Unlike orifice meter, in a nozzle meter the boundary edges are rounded and extended downstream for a short distance. The boundary profile being made to conform to the jet of water coming out of the opening. Thus, there is no contraction of the jet beyond the nozzle opening and therefore coefficient of contraction is unity. Owing to the change in area of flow due to the insertion of the nozzle, there will be difference of pressure between the one point just upstream of the nozzle, there will be difference of pressure between the one point just upstream of the nozzle and one immediately downstream of it. The difference of pressure can be measured by a differential manometer.

To get the relationship for the discharge through a nozzle meter, apply Bernoulli's equation to a point just upstream of nozzle and to a point just downstream of nozzle.

$$\frac{v_1^2}{2g} + \frac{p_1}{\gamma} + z_1 = \frac{v_2^2}{2g} + \frac{p_1}{\gamma} + z_2 \quad (1)$$

The continuity equation is written as, $A_1V_1=A_2V_2$ (2)

Like in case of venturimeter, the above two equations can be solved to give the relation for discharge through a nozzlemeter as

$$Q_{act.} = \frac{C_d A_2}{\sqrt{1 - \left(\frac{d_2}{d_1}\right)^4}} \times \sqrt{2g \left(\frac{p_1 - p_2}{\gamma}\right)} \quad (3)$$

C_d is obtained by calibration of the meter.

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
Procedure

1. Note the dimensions of the cross section of collecting tank and nozzle opening.
2. Keep the drain cock of the collection tank open and open the valve on the pipeline partly to allow some flow through the nozzle meter. Open the flushing cocks of the manometer gradually and simultaneously to remove entrapped air bubbles in the manometer leads and after flushing close the flushing cocks.
3. Note down the reading of the manometer and temperature of the water.
4. Close the drain cock of the collecting tank. The water level in the tank rises. When the level has reached a convenient full scale reading on the scale attached to the gauge glass of the tank, start the stop watch and note the time taken for the known rise of water level in the collecting tank. This completes one set of readings.
5. Vary the discharge through the nozzle meter by manipulating the valve on the pipe line. Allow the condition to become steady. Note the manometer readings and the time for the water level in the collecting tank to rise by a known height after lowering the water level in the collecting tank to a convenient level if necessary by opening the drain cock.
6. Repeat the experiment with various discharges to obtain at least ten sets of readings.
7. Close the valve on the upstream side of the nozzle meter completely, thus cutting out the flow. Close the cocks of the nozzle meter.

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Calculations:

1. Calculate the area of cross section of the collecting tank and area A_2 of the nozzle meter.
2. From the area of cross section of the collecting tank and observed time for water level in the tank to rise by a known height, calculate the discharge, Q for each of the readings.
3. Calculate the values of coefficient of discharge C_d for each set of observations.
4. Calculate Reynolds number for each discharge as $Re = \frac{V_2 d_2}{\nu}$

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6. HYDRAULICS COEFFICIENTS OF ORIFICE METER

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Date / दिनांक: _____

Objective

To determine the hydraulic coefficients of given orifice like Coefficient of Discharge (C_d), Coefficient of Velocity (C_v) and Coefficient of Contraction (C_c)

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Apparatus Required

1. A balancing tank fitted with an orifice provided with a gauge glass and a scale.
2. A Collecting tank with gauge glass and scale and a drain valve.
3. Meter scale
4. Stop Watch.

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Theory:

Outlets in the form of openings or orifices in the walls of storage reservoirs are to be provided for letting out water or any other fluid stored for various utilities. Almost in every case it will be necessary to assess the quantity of fluid issuing from these openings. Theoretically it is possible to obtain the discharge if the area (a) of the orifice and the head (h) of the fluid above the centre of the orifice causing the flow are known, because applying the energy principle to a point in the reservoir and another on the jet issuing from the orifice

$$h = \frac{V_{th}^2}{2g} \text{ or } V_{th} = \sqrt{2gh} \text{ where } V \text{ is the velocity of the jet.}$$

Theoretical discharge $Q_{th} = a\sqrt{2gh}$. But as the streamlines cannot take a sharp turn, they converge to a minimum section of the jet called vena- contracta of area a_c at a distance of about half the diameter of the orifice from the plane of the sharp-edged orifice. The ratio $C_c = (a_c/a)$ is called the coefficient of contraction or $a_c = C_c \cdot a$. Because of loss of head due to resistance to flow by the orifice edge and the change in a direction, the actual velocity V_a of the jet will be less than the theoretical velocity.

The ratio $C_v = \frac{V_a}{\sqrt{2gh}}$ coefficient of velocity or $V_a = C_v \sqrt{2gh}$

Because of the contraction of the jet and the loss of head due to resistance to flow by the orifice boundaries, the actual discharge Q_a will be less than the theoretical

$$\begin{aligned} Q_a &= C_D Q_{th} = a_c V_{th} \\ &= C_c \cdot a \cdot C_v \sqrt{2gh} \\ &= C_d \cdot a \cdot \sqrt{2gh} \text{ Where, } C_d = C_c \cdot C_v \end{aligned}$$

If two of the three coefficients are determined directly, the third can be calculated. C_d is determined in the lab by measuring the actual discharge collected. C_c can be determined by the contraction gauge micrometer but this determination is usually inaccurate. Hence, C_v is calculated from the theoretical velocity, V_{th} and the actual velocity, V_a which is determined from the measured and plotted mean profile of jet as

$V_a = x / \sqrt{2gy}$; x and y being the horizontal and vertical co-ordinates of any point on the mean profile of the jet.

These coordinates of the jet are measured by scale and sliding apparatus (see fig) mounted on the balancing tank vertically above the orifice.

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Procedure

1. Note the dimensions of the collecting tank and the orifice.

2. Open the valve of inflow pipe and adjust the valve such that head in a tank over orifice remain constant while water is flowing from it.
3. Note the head over the orifice and measure x and y coordinates of jet at particular location.
4. Collect water in measuring tank for appropriate time and record the rise of water level in the tank.
5. Repeat steps 2-4.
6. Close the supply valve and drain out the tank and measuring tank and adjust the head over the orifice and repeat the experiment.

Formula used:

Actual discharge - $Q_a = (A*Y)/t$

Theoretical discharge $Q_{th} = a\sqrt{2gh}$

Coefficient of Discharge $C_d = (Q_a / Q_{th})$

Coefficient of Velocity $C_v = \sqrt{x^2/4yh}$

Coefficient of Contraction $C_c = (C_d/C_v)$

Where A - area of measuring tank

a - area of orifice

d - diameter of orifice

H- height of liquid above the centre of the orifice in the tank in which orifice is fitted

t- time to collect water

Y- Rise of water level in piezometer

x- horizontal coordinate of the jet at the measuring point

y- vertical coordinate of the jet at the measuring point

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Calculations

1. Calculate the area of the collecting tank and the area of the orifice.
2. From the area C.S, of the collecting tank and the observed time for the water level in the collecting tank to rise by a known ht. calculate the actual discharge Q_a for each set of readings.
3. Correct the readings of head of water taken on the scale attached to the gauge glass for difference in elevation between the centre of the orifice and the zero of scale. From the corrected head h calculate the theoretical velocity,
$$V_{th} = \sqrt{2gh}$$
4. Calculate the theoretical discharge from the theoretical velocity and orifice area.
5. Calculate $C_d = Q_a/Q_{th}$.
6. Plot the x and y co-ordinates of the jet, referred to the origin at the cr. of the orifice, and draw the mean profile of the jet, taking a point on each of the mean profiles calculate the actual velocities and hence C_v .
7. Calculate $C_c = C_d/C_v$ for each set of readings.

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Graphical Representation of Results:

Draw the curves between

- (1) C_d , C_v , C_c and with head and discharge and
- (2) h with Q_a .

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Questions

1. Where do orifices find application in Engineering works? Explain with a few examples.
2. Explain with a neat sketch how the streamlines passing through an orifice contract to a minimum section at 'Vena Contracta'.
3. Does C_v have greater value for a sharp-edged orifice or a flat and broad edged orifice? Why?
4. Explain with sketches another method of measuring C_v .
5. What will be C_c in round edged orifice?

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Date/ दिनांक.....

Objective: To compute the friction factor for commercial pipe.

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Apparatus required

1. Overhead tank supplying water in which the level is maintained constant.
2. Pipeline with bends valves etc, each fitted with inverted U-tube manometer to measure head loss.
3. A collecting tank and stop watch, to measure the rate of flow.

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Theory:

Loss of head (energy per unit wt. of fluid, expressed in meters of fluid) occurs in any fluid flow. The loss is caused by (i) fluid friction in pipe or conduit (ii) changes in velocity or direction of flow. Losses of these two types are termed as major losses and minor losses respectively. Fluid pipe friction loss is a continuous loss of head, h_f assumed to occur at a uniform rate along the pipe as long as the size and roughness of the pipe remains constant and is commonly referred

to as the loss of head due to friction. For turbulent flow in pipes, the major loss (h_f) is normally taken as

$$h_f = \frac{f l V^2}{2 g d}$$

Where,

f is called the friction factor which is dependent upon the material, surface finish, age etc of the pipe and the Reynolds number of the flow.

l is the length of pipe line in which the loss occurs d is the diameter of pipe line

V is the velocity of flow in the pipeline

Minor losses are losses due to (i) change in cross-section i.e. Expansion loss and Contraction loss (ii) obstruction to flow such as gates, valves etc. (iii) bends, curves etc. Minor losses are generally expressed as a constant time the velocity head, $kV^2/2g$; k differs from fitting to fitting.

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Procedure:

1. Note the diameter (d) and length (l) of the pipe between manometer tapings.
2. Open the valve to allow water to flow through the pipe.
3. Connect all manometers on the pipe line to the flow opening the connecting valves slowly and simultaneously (sudden opening of the valves will cause the water to rush into the manometer limbs and drive the mercury out in case mercury manometers are used).
4. Note down the difference in levels in the columns of the manometer.
5. Change the valve opening and repeat the above steps.
6. Note down the temperature of water.

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Calculations:

Head loss over a length of pipe or in any fitting is given by the manometer connected. From the above equation the friction factor f is calculated.

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Date/ दिनांक.....

Objective: To study various types of flow and to classify them based on Reynolds number.

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Apparatus required:

1. Reynolds apparatus consisting of a balancing tank from which a transparent pipe takes off at the bottom. A regulating valve is fitted at the end of the pipe.
2. A pot of indicator (coloured) liquid with a cock on a thin tube ending into a fine nozzle in the transparent pipe.
3. A collecting tank (measuring jar).
4. Stop watch.

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Theory:

Flow of real fluids in any system offers resistance to flow due to viscosity of the fluid. When the velocity of flow is low, the fluid particles move along smooth lines. Viscous forces are the dominant and curb any tendency in the fluid particle to deviate from their smooth path. If the velocity of the fluid increases after some limit the flow path of the fluid particles becomes wavy as particles tend to deviate from their smooth path. But still the paths are clear. When the velocity is further increased, the fluid particles mix up thoroughly and the flow paths are highly and unidentifiable.

Osborne Reynold studied the nature of flow in pipe line and classified them into various pipes based on a number called 'Reynolds number'.

$$Re = \frac{\rho V D}{\mu} = \frac{V D}{\nu}$$

Where, V is mean velocity of flow, D is diameter of pipe, ρ is the density of fluid flowing, μ is the dynamic viscosity of the fluid and ν is the kinematic viscosity of fluid. When Re is less

than 2000, the flow is laminar. Above this value and up to 4000, the flow is transitional i.e. changing from laminar to turbulent. When Re greater than 4000, the flow is turbulent.

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Procedure:

1. Open the supply valve and admit flow into the balancing tank. When the water level rises adequately, regulate the inflow and outflow so as to keep very low velocity of fluid flowing in the transparent pipe. Keep the head in the balancing tank steady.
2. Admit the indicator fluid (colour liquid) in to the transparent pipe by opening the valve. Observe the flow line of the colour liquid and note the flow pattern.
3. Collect the flow in to the collecting tank (measuring jar) and note the time taken for the volume of water collected.
4. Increase the flow of fluid so that the velocity of flow increases. Again, observe the nature of flow of colour liquid. Observe and note the nature of colour streak.
5. Repeat the above steps to take at least two observations of straight, wavy and zig-zag (mixing) patterns of colour liquid in the transparent tube along with the main flow.
6. Note the diameter of the transparent tube and the temperature of the water.

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Calculations:

1. Calculate the discharge Q , from the volume of water collected in the measuring tank and the time of collection.
2. Calculate the velocity V , of flow through the transparent tube by dividing the discharge by the cross-sectional area of the tube.
3. Calculate the Reynolds number.

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Result:

1. Record the observed nature of the colour streak and the corresponding value of Reynolds number.
2. Write the classification of the flow.

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Date/ दिनांक.....

Objective: To calibrate the water meter.

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Apparatus required:

1. Water meter fitted in a pipe line.
2. Collecting tank, with drain cock, fitted with gauge glass and scale.
3. Stop watch and scale

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Theory:

The water meters are used for knowing the volume of water passed through it over any given period of time. Water meter is of two types.

1. Positive displacement type
2. Inferential type

The positive displacement type has pistons which are displaced by the flowing fluid and a counting mechanism records the number of displacements, The primary parts of an inferential meter are vanes or blades fixed to a shaft or spindle. The shaft rotates due to reaction of the flow on the blades.

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Procedure:

1. Note the dimensions (Cross section) of the collecting tank.
2. Set the inlet regulating valve to an opening thus passing certain discharge through the water meter.
3. Start the stop watch when the water meter pointer comes to a full division on the dial and stop it when the pointer comes to full division again after some time. Note the water meter readings (initial and final) and the time for the indicated flow.
4. Close the outlet valve of the collecting tank and note down the time required for a fixed rise of water level in the collecting tank.
5. Repeat the measurements for various different flow rates by regulating the control valve provided at the end of the pipe line.

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Calculation

Calculate the indicated volume of water flow through the water meter and then the discharge by dividing the volume by the time taken for the indicated flow.

$$Q_{ind.} = V/t.$$

Calculate the actual discharge from the volume of water collected in the collecting tank and time taken for the same. $Q_{act.} = [(L \times B) \times h] / t.$

Calculate the difference between the indicated discharge and the actual discharge as error. -ve of this error is correction to the meter reading.

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Graphical representation

Draw on a simple graph the curve connecting Indicate discharge and actual discharge.

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THIS LAB MANUAL IS PREPARED BY

Dr. Rutuja Manohar Chavan
Assistant Professor

DEPARTMENT OF CIVIL ENGINEERING

Maulana Azad National Institute of Technology, Bhopal

(An Institute of National Importance-Under MHRD, Government of India)

BHOPAL – 462003 (M.P.), INDIA

July 2024